

Brain Health Registry MemTrax Test Information

The MemTrax test is a continuous recognition task that shows 25 unique images and 25 repeat images for a total of 50 images. Participants are instructed to respond only when they see a repeat image. After a participant views the 50th image, the system reports their results.

Actual presentation or performance on this test can be affected by device and internet settings, and all data should be checked before being used for any analyses.

Please note the following when checking the data:

1. Do not use data from an individual performance if there are more than 50 total images.
2. Data collected between 7/4/2018 and 5/31/2019 had only 49 total images (24 unique images and 25 repeat images). Accordingly, the paradigm may have a 2% inconsistency in correct performance metrics compared with tests in which there were 25 unique images. The response times, percentages, and all other values are correct.
3. Data with less than 5 responses or more than 45 total responses should be discarded. These values usually represent individuals who were not responding at all or were responding to all the images.
4. The limits of reasonable response times are 0.3 to 2.3 seconds. Of first-time test takers, 99% are between 0.5 to 1.5 seconds. Repeat test takers may be faster. Rare individuals are slower.
5. For normal and MCI individuals, percent correct should be 70% or more. Of first-time test takers, 99% are better than 75% correct.
6. For clinical populations, MemTrax data is currently being analyzed for the relationship of its metrics to other tests.

